Make the Right Decisions with Reliable Chemical Identification

Basic FT-IR for High-Productivity QA/QC









Long term reliability and tremendous performance in challenging environments enable the Thermo Scientific™ Nicolet™ iS™5 FT-IR spectrometer to serve your QA/QC needs for liquids, solids and gases across polymer, pharmaceutical and many other industrial settings. The Nicolet iS5 spectrometer is a complete tool, offering:

- Full mid-IR spectral range (7800-350 cm⁻¹)
- Tight seals enabling use in humid or dusty locations
- System qualification package for regulatory compliance
- · Lightweight portability and easy maintenance

Fisher Catalog No.	Product Description	Price
14200121	Nicolet iS5 FT-IR spectrometer with OMNIC software, KBr Windows	\$15,500
14200123	Nicolet iS5 FT-IR spectrometer with OMNIC software, ZnSe Windows	\$16,000

Sampling and Software Capabilities



Respond to new analytical demands and improve laboratory productivity with sampling accessories and software that make analysis easy:



- Attenuated Total Reflectance (ATR)
- FT-IR transmission cells and holders
- Thermo Scientific™ OMNIC™ Specta™ mixture analysis software and extensive compound libraries

NEW
Diamond ATR
Accessory



Applications Show You How



Learn how to accomplish more with Nicolet iS5 spectrometer application notes and specialty packages:

Hydrocarbon Levels in Water and Soil

Contaminant Identification

Have questions and want to learn more? Please contact: Contact Information



www.fishersci.com/spectroscopy

©2014 Thermo Fisher Scientific Inc. All rights reserved.
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries

In the United States:

For customer service, call 1-800-766-7000 To fax an order, use 1-800-926-1166 To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437 To fax an order, use 1-800-463-2996 To order online: www.fishersci.ca



A Thermo Fisher Scientific Brand